

# Analysis of plating baths using the Thermo Scientific iCAP 7400 ICP-OES

Patrícia Coelho, Applications Chemist, Thermo Fisher Scientific, Cambridge, UK  
Grégory Lecornet, Applications Chemist, Thermo Fisher Scientific, Courtaboeuf, France

## Key Words

Method of standard addition, plating baths, sulphuric acid levels

## Goal

This application note describes the performance of the Thermo Scientific™ iCAP™ 7400 ICP-OES Duo analysing different elements in different types of plating baths. The duo viewing offers optimal method conditions using axial view for traces and radial view for major elements, and the pre-loaded template allows for rapid and simple method development.

## Introduction

The plating of metal is an ancient technique, used for hundreds of years. It can be defined as the act of covering the surface of objects by depositing a metal on a conductive surface. This is typically carried out by immersing the object in a solution in which the metal ions are moved by an electric field. The applications are rather vast, ranging from the merely decorative, to the enhancement of the physical properties of the material being covered (i.e. to prevent corrosion, reduce friction, alter conductivity, or to improve characteristics such as hardness or durability). It is widely used in the production of jewellery pieces to achieve a silver or gold finish. Due to the ability to cover objects as small as an atom this technique has a high potential to be used in the nanotechnology area.

Manufacturers frequently establish optimum specifications in order to guarantee maximum solution efficiency and uniformity of the plating solution. Plating solutions should be analysed regularly in order to maintain the recommended formulations and to prevent problems related to improper levels of bath constituents and contaminants. A current problem faced by the plating industry are gradual and continuous changes in pH and metal or cyanide content, leading to a significant decrease in efficiency. Some of the most common plating solutions and their main applications are presented in table 1.



Table 1. Plating solutions commonly used and their main applications

Plating Baths	Main applications
Alloy	Hardness improvement
Cadmium	Corrosion Resistance
Chrome	Decorative and industrial
Composite	Physical properties improvement
Golden	Jewellery manufacturing and electronics
Nickel	Decorative and corrosion resistance
Rhodium	Jewellery manufacturing
Silver	Jewellery manufacturing and electronics
Tin	Food processing and electronics
Zinc	Corrosion resistance

There are several different methods that can be employed for the quantitative analysis of plating solutions; these are classified as volumetric, gravimetric or instrumental. Both volumetric and gravimetric methods are simple, accurate, and rapid and can be performed with common laboratory equipment. Nevertheless they only rely on chemical reactions instead of measuring the physical properties associated with the composition of the substance. Additionally, instrumental methods are far quicker and allow for the automatization of the analysis, leading to less mathematical errors and higher reproducibility.

A common instrumental technique is spectroscopic analysis, in particular inductively coupled plasma-optical emission spectroscopy (ICP-OES), which is used in the analysis of major components and trace contaminants in plating solutions. When using this technique there are some aspects that need to be taken into consideration, such as physical interferences (such as viscosity or surface tension) and chemical interferences. They can be overcome easily by sample dilution and accurate matrix matching.

The analysis of sulphuric acid ( $H_2SO_4$ ) levels in plating baths is commonly performed by electroanalytical methods, specifically potentiometry. This is a simple and relatively low cost method but the sensitivity is limited at low concentrations, not allowing for accurate measurements in highly dilute solutions and it can also encounter a number of interferences. Manufacturers that use plating in their processes, are increasingly using ICP-OES instrumentation to quantify  $H_2SO_4$  levels. This method allows for a rapid, sensitive and interference free measurement of sulphur, directly proportional to  $H_2SO_4$ , resulting in high accuracy even at low levels.

## Instrumentation

The Thermo Scientific iCAP 7400 ICP-OES Duo was used for the analysis of a range of plating baths. This is a compact duo view ICP-OES instrument based on the innovative technologies of the iCAP 7000 Series ICP-OES spectrometers. The instrument achieves powerful analyte detection and provides a highly cost effective solution for routine analysis of liquids in laboratories with standard sample throughput requirements. The Thermo Scientific™ Qtegra™ Intelligent Scientific Data Solution™ (ISDS) incorporates several pre-loaded templates (see Figure 1) for common methods, simplifying normal method development and providing an option of immediate analysis.




Figure 1. Plating Baths template selection

## Sample and standard preparation

Elements analysed in each plating bath are presented in table 2. For the Chromium and Chromium & Aluminium baths the quantification of sulphur was performed in order to quantify the amount of  $H_2SO_4$  in the samples.

Table 2. Elements analysed in each plating bath

Type of plating bath	Elements analysed
Chromium and ammonium fluoride (Cr & $NH_4F$ )	Al, Cr, Zr
Nickel (Ni)	B, Ni
Chromium (Cr)	Fe, S
Chromium (Cr)	S
Chromium and aluminium (Cr & Al)	S

Quantification of the samples was made by Method of Standard Addition (MSA) to avoid any possible matrix effect. Standard addition calibration was prepared by spiking samples with traceable 1000mg/L aqueous, single element standards, to the concentrations listed in table 3. All solutions were made to 50 mL with ultra pure de-ionized water.

Table 3. Standard addition concentrations (mg/L)

Type of Plating Bath	Element	Blank	Standard 1	Standard 2	Standard 3
Cr & NH <sub>4</sub> F	Al	0	1	2	-
	Cr	0	20	40	-
	Zr	0	50	100	-
Ni	B	0	5	10	-
	Ni	0	50	100	-
Cr	Fe	0	50	100	-
	S	0	50	100	-
Cr	S	0	50	100	150
Cr & Al	S	0	25	50	75

## Method Development

A LabBook was created in Qtegra ISDS which contained the method parameters and standard concentrations as listed in this note. A standard sample introduction kit was used for the analysis. The instrument was calibrated and the samples analyzed in a single run. The method parameters are shown below in table 4.

Table 4. Method parameters

Parameter	Setting
Pump tubing	Sample Tygon® 1.016 mm Drain Tygon® 1.524 mm
Pump speed	50 rpm
Nebulizer	Glass concentric
Nebulizer gas flow	0.5 L/min
Spray chamber	Glass cyclonic
Center tube	2 mm
RF Power	1150 W
Coolant gas flow	12 L/min
Auxiliary gas flow	0.5 L/min
Exposure times	Low 10 sec High 5 sec

## Results

The results obtained in the analysis of the different bath samples are shown in table 5. It is assumed that all of the sulphur present in the sample is in the form of H<sub>2</sub>SO<sub>4</sub>, therefore the results obtained for sulphur is multiplied by 3.06 [ $M(H_2SO_4)/M(S)=3.06$ ], in order to calculate the concentration of H<sub>2</sub>SO<sub>4</sub>. Results obtained for all the elements were within the expected range.

Table 5. Results of the analysis of the different plating baths. All concentrations are in mg/L.

Type of plating bath	Elements and wavelengths (nm)	Concentration found	Established concentration/range	Dilution factor*
Cr & NH <sub>4</sub> F	Al 167.079	2.2	<10	100
	Cr 205.560	456.0	500	
	Cr 206.550	454.9		
	Cr 267.716	446.7		
	Zr 274.256	1504	1700	
	Zr 327.305	1486		
	Zr 339.198	1567		
	Zr 343.823	1689		
Ni	B 208.959	6380	<7000	1000
	B 249.678	6270		
	B 249.773	6150		
	Ni 221.647	95470	80000	
	Ni 231.604	94770		
Cr	Fe 259.837	2892	<5000	100
	Fe 259.940	2907		
	Fe 371.994	2823		
	S 180.731 (H <sub>2</sub> SO <sub>4</sub> )	1021.6 (3126)	2500 – 3000 (H <sub>2</sub> SO <sub>4</sub> )	
	S 182.034 (H <sub>2</sub> SO <sub>4</sub> )	1035.0 (3167)		
	S 182.624 (H <sub>2</sub> SO <sub>4</sub> )	983.7 (3010)		
Cr	S 180.731 (H <sub>2</sub> SO <sub>4</sub> )	85.6 (2020.3)	2250 – 2500 (H <sub>2</sub> SO <sub>4</sub> )	10
	S 182.034 (H <sub>2</sub> SO <sub>4</sub> )	79.4 (2429.6)		
	S 182.624 (H <sub>2</sub> SO <sub>4</sub> )	74.2 (2271.7)		
Cr & Al	S 180.731 (H <sub>2</sub> SO <sub>4</sub> )	5.7 (35.1)	<150 (H <sub>2</sub> SO <sub>4</sub> )	2
	S 182.034 (H <sub>2</sub> SO <sub>4</sub> )	7.5 (45.9)		
	S 182.624 (H <sub>2</sub> SO <sub>4</sub> )	7.5 (45.9)		

\* All dilutions were prepared using ultra pure water.

Concentration ranges were established by the manufacturers taking into account the concentrations usually found in these types of plating baths to obtain maximum efficiency. This gives the plating operators the status of the bath simplifying the assessment of how the bath composition evolves overtime.

## Conclusion

The analysis of plating baths, specifically those that require the analysis of sulphuric acid, is rapid and highly sensitive when using the Thermo Scientific iCAP 7400 ICP-OES in conjunction with the Qtegra ISDS. This enables the accurate quantification of metals and sulphuric acid content, while the powerful, easy to use iCAP 7400 ICP-OES allows both experienced and inexperienced users alike to vastly reduce the method development time required for these sample types, resulting in cost effective analyses.

[www.thermoscientific.com](http://www.thermoscientific.com)

© 2013 Thermo Fisher Scientific Inc. All rights reserved. Tygon is a registered trademark of Saint-Gobain. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.



Thermo Electron Manufacturing Ltd  
(Cambridge & Courtaboeuf)  
is ISO Certified.

**Africa-Other** +27 11 570 1840  
**Australia** +61 3 9757 4300  
**Austria** +43 1 333 50 34 0  
**Belgium** +32 53 73 42 41  
**Canada** +1 800 530 8447  
**China** +86 10 8419 3588  
**Denmark** +45 70 23 62 60

**Europe-Other** +43 1 333 50 34 0  
**Finland/Norway/Sweden**  
+46 8 556 468 00  
**France** +33 1 60 92 48 00  
**Germany** +49 6103 408 1014  
**India** +91 22 6742 9434  
**Italy** +39 02 950 591

**Japan** +81 45 453 9100  
**Latin America** +1 561 688 8700  
**Middle East** +43 1 333 50 34 0  
**Netherlands** +31 76 579 55 55  
**New Zealand** +64 9 980 6700  
**Russia/CIS** +43 1 333 50 34 0  
**South Africa** +27 11 570 1840

**Spain** +34 914 845 965  
**Switzerland** +41 61 716 77 00  
**UK** +44 1442 233555  
**USA** +1 800 532 4752

AN43150\_E 01/13C

**Thermo**  
SCIENTIFIC  
Part of Thermo Fisher Scientific