thermo**scientific**

Food contact material and food
packaging applications

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**Dionex ASE 350 Accelerated Solvent Extractor**

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Accelerated solvent extraction

Comparison of Soxhlet and accelerated solvent extraction

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
Schematic of the accelerated solvent extraction technique and a Thermo Scientific™ Dionex™ ASE™ 350 Accelerated Solvent Extractor.

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
Chromatograms of MHE extractions from real samples of transparent plastic film.

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
Trace elemental analysis

Peer-reviewed journals

VOCs detected in croissant packaging and their relative distribution.

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
Chromatogram of spiked paperboard with 12 packaging migrants (c = 0.024–30 mg/kg).

Table of contents

Complete FCM and food-packaging solutions

Accelerated solvent extraction

▼ GC and GC-MS

Residual solvents in packaging

VOCs in food packaging

Determination of phthalates

Migration contaminants in paperboard

Phthalate esters in soft drinks

Chemical migrants in breakfast cereal

GC-HRAM


HPLC and UHPLC

LC-MS

LC-HRAM

Trace elemental analysis

Peer-reviewed journals

TIC of 1000 ng/mL phthalate esters standard in full scan 40–450 amu.

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
Table of contents

Complete FCM and food-packaging solutions

Accelerated solvent extraction

GC and GC-MS

Residual solvents in packaging

VOCs in food packaging

Determination of phthalates

Migration contaminants in paperboard

Phthalate esters in soft drinks

Chemical migrants in breakfast cereal

GC-HRAM


HPLC and UHPLC

LC-MS

LC-HRAM

Trace elemental analysis

Peer-reviewed journals

SIM chromatogram for benzophenone and 4-hydroxybenzophenone spiked at 0.6 mg/kg in breakfast cereal.

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Q Exactive GC Orbitrap GC-MS/MS system

Table of contents

Complete FCM and food-packaging solutions

Accelerated solvent extraction

GC and GC-MS

GC-HRAM

Characterizing unknowns in food packaging

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
Zoomed region showing the six peaks of interest in the electron impact (EI) total ion chromatogram of the packaging sample.

Table of contents
Complete FCM and food-packaging solutions
Accelerated solvent extraction
GC and GC-MS
GC-HRAM
HPLC and UHPLC
Determination of bisphenols
Determination of phthalates
LC-MS
LC-HRAM
Trace elemental analysis
Peer-reviewed journals

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
Vanquish UHPLC system

Table of contents

Complete FCM and food-packaging solutions

Accelerated solvent extraction

GC and GC-MS

GC-HRAM

▼ HPLC and UHPLC

Determination of bisphenols



Determination of phthalates

LC-MS

LC-HRAM

Trace elemental analysis

Peer-reviewed journals

**BISPENOL A****BISPENOL B****Chemical structures of Bisphenols A and B.**

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Table of contents

Complete FCM and food-packaging solutions

Accelerated solvent extraction

GC and GC-MS

GC-HRAM

▼ HPLC and UHPLC

Determination of bisphenols

Determination of phthalates

LC-MS


LC-HRAM

Trace elemental analysis

Peer-reviewed journals

Peaks: 1. DMP 2. BMEP 3. DEP
 4. BEEP 5. DPP 6. BBP
 7. DIBP 8. DBP 9. BBEP
 10. DAP 11. DCP 12. BMPP
 13. DHP 14. DEHP 15. DNOP
 16. DNP

Chromatograms of phthalates listed in GB/T 21911-2008 using five different columns.

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Table of contents

Complete FCM and food-packaging solutions

Accelerated solvent extraction

GC and GC-MS

GC-HRAM

HPLC and UHPLC

LC-MS

Analysis of plasticizer contaminants

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
Schematic of method.

Table of contents

Complete FCM and food-packaging solutions

Accelerated solvent extraction

GC and GC-MS

GC-HRAM

HPLC and UHPLC

LC-MS

LC-HRAM

Identification of known and unknown food contact materials

Phthalate screening in food packaging

Trace elemental analysis

Peer-reviewed journals

Vanquish UHPLC and Thermo Scientific™ Q Exactive™ Focus Hybrid Quadrupole-Orbitrap™ MS System

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
Overview of the identified four unknown compounds.

Table of contents

Complete FCM and food-packaging solutions

Accelerated solvent extraction

GC and GC-MS

GC-HRAM

HPLC and UHPLC

LC-MS

▼ LC-HRAM

Identification of known and unknown food contact materials

Phthalate screening in food packaging

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List of analyzed phthalate standards. Exact masses of [M+H]⁺ for precursors as well as for a selection of characteristic HCD fragments are displayed.

Compound	Elemental composition	Precursor [M+H] ⁺	Selection of characteristic HCD fragments m/z	Regulation
Di-n-butyl phthalate (DBP)	C ₁₈ H ₂₂ O ₄	279.1591	167.0339; 205.0859; 223.0965	CA Prop 65
Diisobutyl phthalate (DIBP)	C ₁₈ H ₂₂ O ₄	279.1591	167.0339; 205.0859	EU
Benzyl butyl phthalate (BBP)	C ₁₉ H ₂₀ O ₄	313.1434	91.0542; 205.0859	CA Prop 65
Bis(2-ethylhexyl)phthalate (DEHP)	C ₂₄ H ₃₈ O ₄	391.2843	167.0339; 279.1591	CA Prop 65
Di-n-octyl phthalate (DnOP)	C ₂₄ H ₃₈ O ₄	391.2843	167.0339; 261.1485	CA Prop 65
Diocyltere phthalate (DOTP)	C ₂₄ H ₃₈ O ₄	391.2843	167.0339; 261.1485; 279.1591	-
Diisononyl phthalate (DINP)	C ₂₆ H ₄₂ O ₄	419.3156	127.1481; 275.1642; 293.1747	CA Prop 65
Diisodecyl phthalate (DIDP)	C ₂₆ H ₄₀ O ₄	447.3469	141.1630; 289.1790; 291.1955; 307.1904	CA Prop 65

Table of contents

Complete FCM and food-packaging solutions

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Table of contents

Complete FCM and food-packaging solutions

Accelerated solvent extraction

GC and GC-MS

GC-HRAM

HPLC and UHPLC

LC-MS

LC-HRAM

Trace elemental analysis

Peer-reviewed journals

Recommended resources

www.thermofisher.com/foodcontactmaterials

Peer reviewed journals:

LC and LC-MS methods

Title	Authors	Publication
Release of non-intentionally added substances (NIAS) from food contact polycarbonate: Effect of aging	Chiara Bignardi, Antonella Cavazza, Carmen Laganà, Paola Salvadeo, Claudio Corradini	Food Control, Volume 71, January 2017, Pages 329–335
UHPLC-high-resolution mass spectrometry determination of bisphenol A and plastic additives released by polycarbonate tableware: influence of aging and surface damage	Chiara Bignardi, Antonella Cavazza, Carmen Laganà, Paola Salvadeo, Claudio Corradini	Analytical and Bioanalytical Chemistry, October 2015, Volume 407, Issue 26, Pages 7917–7924
Targeted and untargeted data-dependent experiments for characterization of polycarbonate food-contact plastics by ultra high-performance chromatography coupled to quadrupole Orbitrap tandem mass spectrometry	Chiara Bignardi, Antonella Cavazza, Claudio Corradini, Paola Salvadeo	Journal of Chromatography A, Volume 1372, 12 December 2014, Pages 133–144
Liquid chromatography–full scan-high resolution mass spectrometry-based method towards the comprehensive analysis of migration of primary aromatic amines from food packaging	Monica Mattarozzi, Francesca Lambertini, Michele Suman, Maria Careri	Journal of Chromatography A, Volume 1320, 13 December 2013, Pages 96–102

Methods

Title	Authors	Publication
The challenge of identifying non-intentionally added substances from food packaging materials: A review	C. Nerin, P. Alfaro, M. Aznar, C. Domeño	Analytica Chimica Acta, Volume 775, 2 May 2013, Pages 14–24
Rapid qualitative analysis of phthalates added to food and nutraceutical products by direct analysis in real time/Orbitrap mass spectrometry	Randy L. Self, Wen-Hsin Wu	Food Control, Volume 25, Issue 1, May 2012, Pages 13–16
PCDD and PCDF levels in paper with food contact	H. Beck, A. Droß, W. Mathar	Chemosphere, Volume 25, Issues 7–10, October–November 1992, Pages 1533–1538

Find out more at www.thermofisher.com/food-safety

ThermoFisher
SCIENTIFIC